

# All Ceramics Cement Guide



## Zirconia (Lava, Procera)

- No Etching or bonding
- Use Self Adhesive resin cement, resin-modified glass ionomer cement, or traditional crown & bridge cement

## Lithium Disilicate (IPS e.max Press & CAD)

- Etch tooth surface
- Rinse and apply bonding agent
- Silanate crown before seating
- Cement with light -cure or dual- cure cement OR use self-adhesive resin cement

## Leucite (Authentic, IPS Empress)

- Etch tooth surface
- Rinse and apply bonding agent
- Silanate crown before seating
- Cement with light-cure or dual- cure cement

## Feldspathic

- Etch tooth surface
- Rinse and apply bonding agent
- Silanate crown before seating
- Cement with light-cure or dual-cure cement